

SCALEX™ 830

ULTRASONIC SCALER



SCALEX™
830

- 30,000 scaling strokes per second, solid state power unit.
- Streamlined compact size and slim lightweight handpiece.
- Innovative water-cooled circuitry designed for greater stability during operation.
- No power loss at maximum pressure.
- Available selection of a wide variety of high quality inserts to correspond with specific clinical applications.
- Unit for 25K scaling strokes is available.
- Unit is available in 110/120V 50/60Hz, or 220/240V 50Hz operation.

ISO 9001
ISO 13485



GMP

DENTAMERICA®

SCALEX™ 830

ULTRASONIC SCALER

※ULTRASONIC INSERTS

Reorder No.80005
Universal Insert 25K
(External Water)



Reorder No.83005
Universal Insert 30K
(External Water)



Reorder No.80006
Universal Insert 25K
(Internal Water)



Reorder No.83006
Universal Insert 30K
(Internal Water)



Reorder No.80007
Chisel Type Insert 25K
(External Water)



Reorder No.83007
Chisel Type Insert 30K
(External Water)



Reorder No.80008
Chisel Type Insert 25K
(Internal Water)



Reorder No.83008
Chisel Type Insert 30K
(Internal Water)



Reorder No.80011
Endodontic Insert 25K
(Canalsonic Tip)



SPECIFICATIONS :

- INPUT : 115V 60Hz or 230V 50Hz 3A
- OUTPUT : AC 24V 2.4A
- FREQUENCY : 25,000 or 30,000 Hz/sec.
- DIMENSION : W140xD140xH50mm
- WEIGHT : 1000g

**Reorder
No.88060**
Foot control



Unique compact foot control designed to be activated at any angle and accessible from 360° degrees.

DENTAMERICA®

DISTRIBUTOR

Lions Dental Supply
951-276-3225
www.LionsDentalSupply.com

DENTAMERICA®