

# SCALEX™ 830

## ULTRASONIC SCALER



SCALEX™  
830

- 30,000 scaling strokes per second, solid state power unit.
- Streamlined compact size and slim lightweight handpiece.
- Innovative water-cooled circuitry designed for greater stability during operation.
- No power loss at maximum pressure.
- Available selection of a wide variety of high quality inserts to correspond with specific clinical applications.
- Unit for 25K scaling strokes is available.
- Unit is available in 110/120V 50/60Hz, or 220/240V 50Hz operation.

ISO 9001  
ISO 13485



GMP

**DENTAMERICA®**

# SCALEX™ 830

## ULTRASONIC SCALER

### ※ULTRASONIC INSERTS

**Reorder No.80005**  
Universal Insert 25K  
(External Water)



**Reorder No.83005**  
Universal Insert 30K  
(External Water)



**Reorder No.80006**  
Universal Insert 25K  
(Internal Water)



**Reorder No.83006**  
Universal Insert 30K  
(Internal Water)



**Reorder No.80007**  
Chisel Type Insert 25K  
(External Water)



**Reorder No.83007**  
Chisel Type Insert 30K  
(External Water)



**Reorder No.80008**  
Chisel Type Insert 25K  
(Internal Water)



**Reorder No.83008**  
Chisel Type Insert 30K  
(Internal Water)



**Reorder No.80011**  
Endodontic Insert 25K  
(Canalsonic Tip)



### SPECIFICATIONS :

- INPUT : 115V 60Hz or 230V 50Hz 3A
- OUTPUT : AC 24V 2.4A
- FREQUENCY : 25,000 or 30,000 Hz/sec.
- DIMENSION : W140xD140xH50mm
- WEIGHT : 1000g

**Reorder  
No.88060**  
Foot control



Unique compact foot control designed to be activated at any angle and accessible from 360° degrees.

**DENTAMERICA®**

DISTRIBUTOR

**Lions Dental Supply**  
951-276-3225  
[www.LionsDentalSupply.com](http://www.LionsDentalSupply.com)

**DENTAMERICA®**