

# SCALEX™ 850

SELF-CONTAINED ULTRASONIC SCALER



Self-contained solution tank

- Self-contained solution tank allows for convenient portability.
- One touch selection between medicament tank and external water source.
- 25k/30k scaling strokes per second, solid state power unit.
- Innovative water-cooled circuitry designed for greater stability during operation.
- Available selection of a wide variety of high quality inserts to correspond with specific clinical applications.
- Unit is available in 100/110V 50/60Hz or 220/240V 50Hz operation.

ISO 9001  
ISO 13485



GMP

**DENTAMERICA®**

P.O. Box 3200  
Industry, CA 91744 U.S.A.

Phone: (626) 912-1388  
Fax: (626) 913-0510  
[www.dentamerica.com](http://www.dentamerica.com)

# SCALEX™ 850

## ULTRASONIC SCALER

### ※ULTRASONIC INSERTS

**Reorder No.80005**

Universal Insert 25K  
(External Water)

**Reorder No.83005**

Universal Insert 30K  
(External Water)

**Reorder No.80006**

Universal Insert 25K  
(Internal Water)

**Reorder No.83006**

Universal Insert 30K  
(Internal Water)

**Reorder No.80007**

Chisel Type Insert 25K  
(External Water)

**Reorder No.83007**

Chisel Type Insert 30K  
(External Water)

**Reorder No.80008**

Chisel Type Insert 25K  
(Internal Water)

**Reorder No.83008**

Chisel Type Insert 30K  
(Internal Water)

**Reorder No.80011**

Endodontic Insert 25K  
(Canalsonic Tip)

**Reorder No.88060****Foot control**

Unique compact foot control  
designed to be activated at  
any angle and accessible  
from 360° degrees.

**SPECIFICATIONS :**

- \* INPUT : 115V 60Hz or 230V 50Hz 3A
- \* OUTPUT : AC 24V 2.4A
- \* FREQUENCY : 25,000 or 30,000 Hz/sec.
- \* DIMENSION : W140xD140xH50mm
- \* WEIGHT : 1000g

**DENTAMERICA®****Reorder No.85003****Self-contained solution tank**

Easily change the solutions according to different  
applications including periodontics, endodontics,  
and prophylaxis.

DISTRIBUTOR

**Lions Dental Supply****951-276-3225**[www.LionsDentalSupply.com](http://www.LionsDentalSupply.com)**DENTAMERICA®**