

# SCALEX™ 800

## ULTRASONIC SCALER



- 25,000 scaling strokes per second, solid state power unit.
- Streamlined compact size and slim light weight handpiece.
- Convenient foot control.
- No power loss at maximum pressure.
- Unit is available in 115V 60Hz, or 230V 50Hz operation.
- Dimensions W220 x D190 x H80 mm, weight 3.2kg.

ISO 9001  
ISO 13485



GMP

**DENTAMERICA®**

# SCALEX™ 800

## ULTRASONIC SCALER

※ULTRASONIC INSERTS



<b>Reorder No.80005</b>	<b>Reorder No.80006</b>	<b>Reorder No.80007</b>	<b>Reorder No.80008</b>	<b>Reorder No.80011</b>
Universal Insert 25K (External Water)	Universal Insert 25K (Internal Water)	Chisel Type Insert 25K (External Water)	Chisel Type Insert 25K (Internal Water)	Endodontic Insert 25K (Canalsonic Tip)

### SPECIFICATIONS :

- INPUT : 115V 60Hz, or 230V 50Hz
- OUTPUT : DC 24V
- FREQUENCY : 25,000 scaling stroke per second
- DIMENSION : W220HD190HH80mm
- WEIGHT : 3.2kg

**Reorder  
No.80170**  
Variable Speed  
Foot Control



**DENTAMERICA®**

DISTRIBUTOR

**Lions Dental Supply**  
**951-276-3225**  
[www.LionsDentalSupply.com](http://www.LionsDentalSupply.com)

**DENTAMERICA®**