

SCALEX™ 880 PLUS

PIEZO ULTRASONIC SCALER



- State-of-the-art piezoelectric technology is self-tuning and load-adjusting for smooth and effortless operation.
- Special TURBO function helps remove especially hard calculi by increasing power relative to the hardness.
- Streamlined compact size and slim lightweight handpiece.
- Innovative water-cooled circuitry designed for greater stability during operation.
- Extremely gentle on the patient.
- Available selection of a wide variety of high quality tips to correspond with specific clinical applications.
- Unit is available for 100-240V 50/60Hz operation.
- Three tips included.

ISO 9001
ISO 13485



DENTAMERICA®

SCALEX™ 880 PLUS

SCALEX™ 880 PLUS

PIEZO ULTRASONIC SCALER

※PIEZO SCALER TIPS



Reorder No. 88006
Perio



Reorder No. 88004
Universal B



Reorder No. 88005
Chisel



Reorder No. 88003
Universal A



Reorder No. 88002
Endo Chuck

SPECIFICATIONS:

- INPUT: 100V-240V 50/60Hz 0.8A
- OUTPUT: DC 24V 1.5A
- FREQUENCY: 27,000 to 32,000 Hz/sec
- DIMENSIONS: W140xD140xH50mm
- WEIGHT: 1,000g

DENTAMERICA®



Reorder No. 88015
Handpiece

Autoclavable handpiece with piezo ceramic crystals and titanium tip socket structure.



Reorder No. 88060
Foot Control

Unique compact foot control designed to be activated at any angle and accessible from 360° degrees.



Reorder No. 88009
Torque-limiting
Tip Wrench

Special autoclavable wrench prevents damage to scaler tip when tightening.

DISTRIBUTOR

Lions Dental Supply
951-276-3225
www.LionsDentalSupply.com

DENTAMERICA®