

# SCALEX™ 890

SELF-CONTAINED PIEZO ULTRASONIC SCALER



Self-contained solution tank

SCALEX™ 890

- Self-contained solution tank allows for convenient portability.
- One touch selection between medicament tank and external water source.
- State of the art piezoelectric self-tuning and load-adjusting technology for smooth and effortless operation.
- Innovative water-cooled circuitry designed for greater stability during operation.
- Available selection of a wide variety of high quality tips to correspond with specific clinical applications.
- Unit is available in 100/110V 50/60Hz or 220/240V 50Hz operation.

ISO 9001  
ISO 13485



**DENTAMERICA®**

# SCALEX™ 890

## SELF-CONTAINED PIEZO ULTRASONIC SCALER

### ※PIEZO SCALER TIPS



Reorder No.88002  
Endo Chuck



Reorder No.88003  
Universal A



Reorder No.88004  
Universal B



Reorder No.88005  
Chisel



Reorder No.88006  
Perio



Reorder No.89003  
Self-contained solution tank

Easily change the solutions according to different applications including periodontics, endodontics, and prophylaxis.



Reorder No.88015  
Handpiece

Autoclavable handpiece with piezo ceramic crystals and titanium tip socket structure.



Reorder No.88060  
Foot control

Unique compact foot control designed to be activated at any angle and accessible from 360 degrees.

### SPECIFICATIONS:

- \* FREQUENCY: 27,000 to 32,000 Hz/sec.
- \* DIMENSIONS: D260xW140xH210 mm
- \* VOLUME: 1,000 ml
- \* WEIGHT: 3.3 kg
- \* INPUT: 115V/60Hz or 230V/50Hz
- \* OUTPUT: AC 24V 3A

**DENTAMERICA®**

DISTRIBUTOR

**Lions Dental Supply**  
**951-276-3225**  
**[www.LionsDentalSupply.com](http://www.LionsDentalSupply.com)**

**DENTAMERICA®**