



# ROTEX™ 782

## COMPACT ELECTRIC HANDPIECE UNIT

- Ball bearing high torque electric motor.
- 1,000 to 35,000 R.P.M. variable speed.
- Double insulation heavy duty construction.
- Automatic cooling fan for continuous high speed performance.
- Forward/Reverse switch.
- Standard package includes Variable Speed Foot Control, Handpiece and Handpiece Rest.
- Circuit overload protection.
- Application for grinding, polishing, milling of precious metals, ceramics, acrylics and steel.
- Unit is available for 100/110V 50/60Hz or 220/240V 50Hz operation.

**Lions Dental Supply**  
951-276-3225

[www.LionsDentalSupply.com](http://www.LionsDentalSupply.com)

ISO 9001  
ISO 13488



ROTEX™ 782