

ROTEX™ 782E

COMPACT ELECTRIC E TYPE HANDPIECE UNIT



ISO 9001
ISO 13485



- Ball bearing high torque electric motor.
- 1,000 to 35,000 R.P.M. variable speed.
- Double insulation heavy duty construction.
- Automatic cooling fan for continuous high speed performance.
- Forward/Reverse switch.
- Application for grinding, polishing, milling of precious metals, ceramics, acrylics and steel.
- Standard package includes Variable Speed Foot Control, Handpiece and Handpiece Rest. Circuit overload protection.
- Unit is available for 100/110V 50/60Hz or 220/240V 50Hz operation.

DENTAMERICA®

ROTEX™
782E

ROTEX™ 782E

COMPACT ELECTRIC E TYPE HANDPIECE UNIT

Reorder No. 7822

Straight Attachment



Reorder No. 7824

Contra Attachment

Reorder No. 783

E-Type Motor Handpiece



Reorder No. 782208

Variable Speed Foot Control



SPECIFICATIONS :

- INPUT : 115V/230V 50/60Hz
- OUTPUT : DC 2V-39V
- SPEED : 1,000-35,000 rpm
- DIMENSIONS : L120xW190xH80mm
- WEIGHT : 2.0kg

DENTAMERICA®

DISTRIBUTOR

Lions Dental Supply
951-276-3225

www.LionsDentalSupply.com

DENTAMERICA®