

JDS Dental Products

Vortex Series Wet Ring Vacuum Pumps



Single							
MoPel Number	Total HP	User Capacity	Voltage	Line Frequency	Total Amps	Size (H x W x D)	Shipping Weight
JVBS10	1.0	2	115/208-230	60 Hz	16/8	16 x 12 x 12	58 lbs
JVBS10R	1.0	2	115/208-230	60 Hz	16/8	16 x 12 x 12	60 lbs
JVBS10S	1.0	2	115/208-230	60 Hz	16/8	16 x 12 x 12	60 lbs
JVBS10RS	1.0	2	115/208-230	60 Hz	16/8	16 x 12 x 12	62 lbs
JVBS15	1.5	3	115/208-230	60 Hz	22/11	17 x 12 x 12	62 lbs
JVBS15R	1.5	3	115/208-230	60 Hz	22/11	17 x 12 x 12	64 lbs
JVBS15S	1.5	3	115/208-230	60 Hz	22/11	17 x 12 x 12	64 lbs
JVBS15RS	1.5	3	115/208-230	60 Hz	22/11	17 x 12 x 12	66 lbs
JVBS20	2.0	4	208-230	60 Hz	13.4	19 x 12 x 12	75 lbs
JVBS20R	2.0	4	208-230	60 Hz	13.4	19 x 12 x 12	77 lbs
JVBS20S	2.0	4	208-230	60 Hz	13.4	19 x 12 x 12	77 lbs
JVBS20RS	2.0	4	208-230	60 Hz	13.4	19 x 12 x 12	79 lbs

Twin							
MoPel Number	Total HP	User Capacity	Voltage	Line Frequency	Total Amps	Size (H x W x D)	Shipping Weight
JVBTW20S	2.0	4	208-230	60 Hz	8 (2x)	23 x 28 x 17	132 lbs
JVBTW20RS	2.0	4	208-230	60 Hz	8 (2x)	23 x 28 x 17	135 lbs
JVBTW30S	3.0	6	208-230	60 Hz	11 (2x)	26 x 28 x 17	150 lbs
JVBTW30RS	3.0	6	208-230	60 Hz	11 (2x)	26 x 28 x 17	154 lbs
JVBTW40R	4.0	8	208-230	60 Hz	13.4 (2x)	30 x 28 x 17	160 lbs
JVBTW40RS	4.0	8	208-230	60 Hz	13.4 (2x)	30 x 28 x 17	165 lbs

- R Penotes Water Recycler (0.19 GPM)
- S Penotes Air/Water Separator
- Twin anP Triple moPels can be custom maPe to orPer.

Showcasing our most PefineP line of wet-ring vacuums, JDS presents the Vortex Series Wet Ring vacuum pumps. With its Pebut in 2004, the Vortex Series has been rePesigneP for maximum performance anP quality, proviPing years of PepenPable service.

EquipeP with an impressive 3 year warranty, along with a minimal maintenance level, built to last in continuous operation, chemical resistant marine bronze construction, anP rePuceP noise level, all for your comfort in your office.

Key features

- EquipeP with motors from Bluffton Motor Works, reliability anP trouble free operation comes first
- ETL listeP, UL #60601-I Compliant: StanParP for safety mePical electrical equipment
- CertifieP to CSA C22.2 #601.1: MePical electrical equipment
- Water recycler saves 70% on water consumption
- Dual voltage for 1 HP anP 1.5 HP moPels, single voltage for 2 HP moPels



To Order Call
Lions Dental Supply
951-276-3225