

The Dynamic Duo



We have **OPTIMIZED** two of our most popular units into **ONE COMPLETE System!!**

Our V-35 MicroMotor – The V35 electric handpiece system features a High Quality, multi-speed control console teamed with our Dependable Gray Handpiece. A heavy-duty, 35,000 RPM, lab-only handpiece system for use with HP-sized burs that offers high durability, high torque and extremely smooth operation. NO Foot Control Required!!!!



AND

Our ThermaKnife™ Thermal Tray Trimming Knife - A time-saving knife with a hot-tipped blade that promotes fast and easy trimming of ALL vacuum formed materials. Adjustable temperature control for precise heat adjustment providing the operator with immediate heat for fast cutting and trimming. Ultra-fast heating circuit gets tip to ideal temperature in seconds, saving time. Lightweight pen with sharp, permanent tip promotes improved comfort and control. We have even added an EXTRA OUTLET for an additional Pen. Unit Comes complete with 1 Knife edged Tip Pen (additional Pens available – please see current Catalog.)



Ordering Information:

#38760 Buffalo Dynamic Duo V-35/Therma-Knife Combo Unit Complete

To Order Call:

Lions Dental Supply
951-276-3225
www.LionsDentalSupply.com

To Order Call
Lions Dental Supply
951-276-3225
www.LionsDentalSupply.com