

# MARATHON



knee control style  
Mounting bracket for knee control

## Multi 600

- Multi switching function.
  - Easy to change Desktop style into **knee control style** with a mounting bracket
- Speed: 1,000 ~ Max 50,000 rpm.
- Low noise, Low vibration.
- Excellent Durability.
- Ceramic Ball Bearing.
- Ergonomic Design.
- Self-Diagnosis Function.
- Error Code Display Function.
- Auto-Cruise Function.

### Control unit Model. **SDE-S60/L60**

▶ Catalog No. #6000 ▶ Voltage 100~120V/200~240V ▶ Cycle 50/60Hz  
 ▶ Electric Power 100VA ▶ R.P.M. Handpiece Max ▶ Weight 3.1Kg  
 ▶ Dimension 122(W)mm × 268(D)mm × 237(H)mm

Desktop style



### Motor handpiece Model. **BM50S1**

▶ Catalog No. #6102 ▶ nonload 1,000~50,000rpm  
 ▶ Max Torque 7.8Ncm ▶ Electric Power 230W ▶ Weight 238g  
 ▶ Dimension 159.8(L)mm × 26.2(D)mm

## Multi 600 STANDARD

- Vertical Design for desktop
  - Foot or Hand Operated-
- Speed: 1,000 ~ Max 50,000 rpm.

### Control unit Model. **SDE-S60s**

▶ Catalog No. #6001 ▶ Voltage 100~120V/200~240V ▶ Cycle 50/60Hz  
 ▶ Electric Power 100VA ▶ R.P.M. Handpiece Max ▶ Weight 3.1Kg  
 ▶ Dimension 122(W)mm × 268(D)mm × 237(H)mm

### Motor Handpiece Model. **BM50S1**

▶ Catalog No. #6102 ▶ nonload 1,000~50,000rpm  
 ▶ Max Torque 7.8Ncm ▶ Electric Power 230W  
 ▶ Weight 238g ▶ Dimension 159.8(L)mm × 26.2(D)mm



## Multi 600 DUO

- Dual control system both **brushless** and **carbon brush** motor.
- Multi switching function.
  - Easy to change Desktop style into knee control style with a mounting bracket
- Speed: 1,000 ~ Max 40,000 rpm.

### Control unit Model. **SDE-S60/L60D**

▶ Catalog No. #6002 ▶ Voltage 100~120V/200~240V ▶ Cycle 50/60Hz ▶ Electric Power 60VA  
 ▶ R.P.M. Handpiece Max ▶ Weight 2.9Kg ▶ Dimension 122(W)mm × 268(D)mm × 228(H)mm

### Motor Handpiece Model. **BM40S1**

▶ Catalog No. #6212 ▶ nonload 1,000~40,000rpm ▶ Max Torque 8.5Ncm ▶ Electric Power 230W  
 ▶ Weight 238g ▶ Dimension 159.8(L)mm × 26.2(D)mm

