



2005 Orthodontic Patient Chair




Panel Mount 4062 for two handpiece manual
 SP400FB Ortho cabinet with keyboard drawer in fixed base
 FPMA103G Panel mount arm (grey) (monitor not included)
 D5000 Ultra doctor's stool in optional Papyrus Ultraleather
 2005 shown in optional Papyrus Ultraleather Plush upholstery



Panel Mount 4062 for two handpiece manual
 SP200CA Ortho standard mobile cabinet on caster base
 D5000 Ultra doctor's stool in optional Papyrus Ultraleather
 2005 shown in optional Papyrus Ultraleather Plush upholstery



2005 Orthodontic Chair with
 Electro-mechanical backrest



Estimated shipping dimensions and weights
 Chair - 41L x 31W x 28H 240 lbs. plus
 Upholstery - 47L x 11W x 23H 30 lbs.

STANDARD FEATURES	PART NO.	MODEL	2000064T	2000-064S	2000-064SA	2000-064
<i>Suggested Retail Price</i> 2005 Plush with trim Orthodontic patient chair (available in ultraleather only)			\$3,607			
2005 Slings Plush with trim Orthodontic patient chair (available in ultraleather only)				\$3,607		
2005 Slings asepsis Orthodontic patient chair					\$3,171	
2005 Asepsis Orthodontic patient chair						\$3,069
Foot Switches with automatic exit (return), mounted on pedestal			❖	❖	❖	❖
Electromechanical backrest			❖	❖	❖	❖
All-metal construction with HiGlos powdercoat inside and out			❖	❖	❖	❖
Pedestal base with 5" of manual height adjustment			❖	❖	❖	❖
Double-articulating headrest <i>adjustable for height and wheelchair positioning</i>			❖	❖	❖	❖
Synchronized seat and back tilt			❖	❖	❖	❖
Quick change upholstery (Plush style done in Ultraleather only)			❖	❖	❖	❖
Three year warranty			❖	❖	❖	❖
FACTORY OPTIONS	PART NO.	PRICE				
Accessory power cord	2000-238	\$144	❖	❖	❖	❖
Chair mount adapter	2000-226	\$250	❖	❖	❖	❖
Stationary arms	2000-243	\$200	❖	❖	❖	❖
Moveable arms (arms swing up)	2000-244	\$355	❖	❖	❖	❖
Arm covers for Asepsis 2005 with added arm feature	7000-096-1	\$95/set				❖
Touchpad with 5 programmable presets and two automatic exit (return) must also purchase Preset circuit board (see below)						
Chair mounted	2000-231	\$374	❖	❖	❖	❖
Preset circuit board (one preset will work with the pedestal switches, must purchase foot control or touchpad to work more than one preset)	7003-150-2A	\$300	❖	❖	❖	❖
Foot Control - 7 function	7000-060-2	\$350	❖	❖	❖	❖
Non-articulating headrest (can be purchased or in lieu of articulating)	7000-028	\$429	❖	❖	❖	❖
Ponytail headrest (\$100 if purchased in lieu of articulating)	7000-609	\$529	❖	❖	❖	❖
Replacement Double articulating headrest with upholstery	7000-021	\$429	❖	❖	❖	❖
Replacement upholstery, Plush with trim (available in ultraleather only)	7005-703	\$1,438	❖			
Replacement upholstery, Sling plush with trim (available in ultraleather only)	7005-706	\$1,438		❖		
Replacement upholstery, Sling asepsis	7005-709	\$1,000			❖	
Replacement upholstery, asepsis (slide on back)	7005-700	\$900				❖
Replacement headrest upholstery	7005-041	\$80	❖	❖	❖	❖
Ultraleather upgrade (not plush styling)	1005	\$404	❖	❖	❖	❖
Florentine Burl Trim	2000-229	\$215/set	❖	❖	❖	❖
Extended Three Year Warranty for each chair	1008	\$400	❖	❖	❖	❖