

OSIII Oral Surgery Patient Chair



Shown with Eight way I. V. arm option

Shown in Duration Taupe Upholstery



Standard I. V. arm pillow
Interchangeable Left and Right




Eight way articulating I. V. arm
option
Interchangeable Left and Right



Chinner headrest option

STANDARD FEATURES	PART NO.	MODEL	2000-080-2
<i>Suggested Retail Price</i> OSIII Oral surgery electromechanical patient chair			\$6,242
Foot Control with two programmable preset positions and automatic exit (return), 7 function			❖
Electromechanical Cantilever system			❖
All-metal construction with HiGlos powdercoat inside and out			❖
Thermal overload protection			❖
Left/Right I. V. arms with one interchangeable pillow			❖
Double-articulating headrest with lock down for surgery feature <i>adjustable for height and wheelchair positioning</i>			❖
Synchronized seat and back tilt			❖
Quick change upholstery			❖
Upholstered arms			❖
Three year warranty			❖
FACTORY OPTIONS	PART NO.	PRICE	
Touchpad with 5 programmable presets and two automatic exit (return)			
Chair mounted	2000-231	\$374	❖
Delivery mounted	2000-234		
Chinner double articulating headrest (can be purchased as upgrade for \$308)	7000-608	\$737	❖
Eight way articulating I. V. arm	2000-220	\$586	❖
Leg and Body strap (color: black only)	2000-223	\$164	❖
Additional I. V. pillow only	7000-116	\$162	❖
Height spacer 3" (factory installed only)	2000-119	\$205	❖
Mobile locking base	7000-371	\$1,251	❖
Chair mount adapter	2000-226	\$250	❖
Replacement Double articulating headrest with upholstery	7000-021	\$429	❖
Replacement upholstery asepsis	7000-535	\$1,238	❖
Replacement headrest upholstery	7005-041	\$80	❖
Ultraleather upgrade (not plush styling)	1005	\$404	❖
Extended Three Year Warranty for each chair	1008	\$400	❖



Estimated shipping dimensions and weights
 Chair - 41L x 31W x 28H 300 lbs. plus
 Upholstery - 47L x 11W x 23H 30 lbs.