

Lil' BEAVER 2.0

MAGNETOSTRICTIVE Ultrasonic Scaling Unit



Accepts Both
25K and 30K
Inserts

Insert not included



**Includes vertical
mount base!**

- Accepts both 25K (long) and 30K (short) ultrasonic inserts - no adapter or switch, unit automatically adjusts to accommodate any style, brand or frequency.
- Convenient front-mounted water control with large easy-to-adjust dial - easy even with gloves.
- Turbo boost, 2 position foot pedal - press lightly for regular dial power, press firmly for TURBO! *Heavy calculus watch out!*
- Digital power adjustment with LED lit power display panel.
- Front mounted main power switch - shutting the unit off when left unattended just got a whole lot easier.
- All aluminum housing for improved durability. No more plastic to stain, crack or break
- Lies flat or mounts vertically in provided vertical mounting stand for optimal counter space saving. Fits almost anywhere!



Complete Kit Includes:
Little Beaver Scaling Unit, Power Cord, 2 Position Turbo Boost Foot Pedal,
1/4" Water Line w/Quick-Disconnect, Vertical Mount Stand, Handpiece
Holder Mount, Autoclavable Handpiece Sheath. 1 year warranty.
SC-001 Retail \$ **799.00** ea

Lions Dental Supply
951-276-3225
www.LionsDentalSupply.com