




Save up to 50% on this Portable X-ray Generator compared to Aribex Nomad Pro

Why MobileX?


- Stay With Your Patient the Whole Time
- Reduce Retakes
- Ideal for Children & Phobic Patients
- Take X-rays At Any Point In The Procedure
- FDA Registered – K180561
- Film, Plate & Sensor Compatibility
- Smaller, lighter, easier to handle
- Low Patient Exposure
- Camera-like Design
- Portable for Ease of Use
- No Installation
- No Walls Required
- Take X-rays In Any Operatory
- DC Unit
- High Image Quality
- Convenient & Mobile
- Replaces Wall-Mounted Units



The portable x-ray solution for new practices!
A great alternative to bulky wall-mounted x-ray generators!
Replaces your older x-ray heads, with lower settings!
Go digital-fast, easy & economically!

Technical Specifications

X-ray Generator:	Type: High Frequency Inverter
X-ray Tube:	Latest Toshiba Technology
	Tube voltage: 70KV DC
	Tube current: 2mA
	Anode Type: Stationary
	Anode Angle: 12.5°
	Focal Spot: 0.4x0.4 mm
X-ray Control:	Exposure time from 0.01 to 2.0sec (0.01 sec step)
Power Requirement:	12.6 VDC 1.5A
Weight:	5.2 lbs including battery, shield
Dimensions:	6.5"(W) x 6.1"(H) x 2.4"(D)
Battery:	Li-Po 11.1V



The compact "Camera-like" generator increases the user's confidence and improves patient perception of x-ray's. MobileX increases dental radiography speed, convenience, efficiency and production. MobileX exposes both the patient and operator to the least amount of radiation. Just one MobileX does the work of numerous conventional x-ray systems, making it extraordinarily economical.

Lions Dental Supply 951-276-3225